



# PURline™ Composite Plates Compatibility Tool



|                                                  | <b>AZ-7283</b><br>0.25 mL<br>96 Well<br>Half-Skirt<br>Clear | <b>AZ-7284</b><br>0.25 mL<br>96 Well<br>Half-skirt<br>White | <b>AZ-7281</b><br>0.15 mL<br>96 Well<br>Full-skirt<br>Clear | <b>AZ-7282</b><br>0.15 mL<br>96 Well<br>Full-skirt<br>White | <b>AZ-7285</b><br>35 µL<br>384 Well<br>Full-skirt<br>Clear | <b>AZ-7286</b><br>35 µL<br>384 Well<br>Full-skirt<br>White |
|--------------------------------------------------|-------------------------------------------------------------|-------------------------------------------------------------|-------------------------------------------------------------|-------------------------------------------------------------|------------------------------------------------------------|------------------------------------------------------------|
| <b>Applied Biosystems (Thermo)</b>               |                                                             |                                                             |                                                             |                                                             |                                                            |                                                            |
| 3100, 3130XL, 3500, 3500XL, 3730, 3730XL         |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| GeneAmp 2700, 2720                               | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| GeneAmp 9600                                     | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| GeneAmp 9700                                     | ✓                                                           | ✓                                                           |                                                             |                                                             | ✓                                                          | ✓                                                          |
| MiniAmp                                          | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| ProFlex                                          | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| Proflex, Multiblock system                       |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| QuantStudio 5, 6, 7, 12 K Real-time PCR System   | ✓                                                           | ✓                                                           |                                                             |                                                             | ✓                                                          | ✓                                                          |
| SimpliAmp                                        | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| Veriti 384-well                                  |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| Veriti 96-well                                   | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| ViiA7, 7900 HT FAST                              |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| <b>Agilent</b>                                   |                                                             |                                                             |                                                             |                                                             |                                                            |                                                            |
| AriaMx qPCR System                               |                                                             |                                                             | ✓                                                           | ✓                                                           |                                                            |                                                            |
| Mx3000P, Mx3005P                                 | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| Robocycler Gradient                              | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                           |                                                            |                                                            |
| Surecycler 8800                                  | ✓                                                           | ✓                                                           |                                                             |                                                             | ✓                                                          | ✓                                                          |
| <b>BioRad</b>                                    |                                                             |                                                             |                                                             |                                                             |                                                            |                                                            |
| C1000, C1000 Touch, S1000                        | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                          | ✓                                                          |
| CFX Connect, CFX96, CFX96 Touch                  |                                                             |                                                             | ✓                                                           | ✓                                                           |                                                            |                                                            |
| CFX384, CFX384 Touch                             |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| CFX96 Touch Deep Well                            | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                           |                                                            |                                                            |
| Chromo4                                          | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                           |                                                            |                                                            |
| DNA Engine Opticon and Opticon 2                 |                                                             |                                                             | ✓                                                           | ✓                                                           |                                                            |                                                            |
| DNA Engine, DNA Engine Tetrad, Tetrad 2          |                                                             |                                                             | ✓                                                           | ✓                                                           | ✓                                                          | ✓                                                          |
| iCycler                                          | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| iQ, iQ-4, iQ-5, myiQ PCR System                  | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| MJ Research (Bio-Rad) PTC200, 220, 221, 225, 240 | ✓                                                           | ✓                                                           |                                                             |                                                             | ✓                                                          | ✓                                                          |
| MyCycler                                         | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| PTC 200                                          |                                                             |                                                             | ✓                                                           | ✓                                                           | ✓                                                          | ✓                                                          |
| T100                                             | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| <b>Corbett</b>                                   |                                                             |                                                             |                                                             |                                                             |                                                            |                                                            |
| Corbett Research (Qiagen) Palm Cycler            | ✓                                                           | ✓                                                           |                                                             |                                                             |                                                            |                                                            |
| Corbett Research (Qiagen) Palm Cycler 384        |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| <b>Eppendorf</b>                                 |                                                             |                                                             |                                                             |                                                             |                                                            |                                                            |
| MasterCycler ep 384 / Pro 384                    |                                                             |                                                             |                                                             |                                                             | ✓                                                          | ✓                                                          |
| Mastercycler ep realplex                         | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                           |                                                            |                                                            |
| MasterCycler ep, ep gradient, Pro, Pro S, nexus  | ✓                                                           | ✓                                                           | ✓                                                           | ✓                                                           |                                                            |                                                            |